

## Pharmacy Fridge - 66 Litres

### Product Images



### Description

- Compliant with all leading professional guidelines...
  - Department of Health Green Book
  - Recommendations in The Royal Pharmaceutical Society of Great Britain's 'Code of Ethics & Standards'
  - National Travel Health Network and Centre - Yellow Book
  - Royal College of Veterinary Surgeons Guidelines
- Fan assisted cooling
- Designed specifically for the safe storage of vaccines and pharmaceuticals
- Available as bench top, under bench, wall fitting or free-standing
- Temperature range +2°C to +8°C
- Second wandering temperature probe (can be placed anywhere within the refrigerator, ideal for pharmacists needing to comply with MHRA regulations)
- IntelliCold® touch screen controller
  - Automatic 24/7 temperature & fridge operation data recording (micro SD card required)
  - Min/max temperature control, display & recording
  - Displays temperature from both probes
  - High/low temperature alarm
  - Door open alarm
  - Mains failure alarm (requires batteries)
  - Alarm battery back-up
- Internal LED light
- Access port & remote alarm contacts
- Fan assisted cooling
- Off cycle auto defrost
- Low GWP environmentally friendly refrigerant
- Self-closing doors with door open alarm
- Available with either a CAM or digital push button lock
- CAM - Supplied with two keys
- Digital Push Button - Mains operated (uses a 4-digit pin)
- Adjustable shelves
- Micro SD card (optional)
- UKAS traceable calibration service (optional)
- Repeat alarm (optional)
- **Five year parts & labour warranty**

For information on recycling and cleaning instructions please refer to [www.pharmacyequipmentdirect.com](http://www.pharmacyequipmentdirect.com) or contact our technical services department on 01258 451338



## Additional Information

Specification	Model No.	5RLDF0219
	Door Type	Solid
	Capacity	66 Litres
	Dimensions (w x d x h)	450 x 510 x 740mm
	Weight	34.0kg
	Positioning	Bench/Wall Mounting
	Shelves	3
	Drawer Capacity	-
	Energy Consumption	0.63 kWh/24h
	Operating Temp (nominal)	5°C
Temperature Range	2 to 8°C	

For information on recycling and cleaning instructions please refer to [www.pharmacyequipmentdirect.com](http://www.pharmacyequipmentdirect.com) or contact our technical services department on 01258 451338

